FOR IN VITRO RESEARCH USE ONLY (RUO) NOT FOR CLINICAL DIAGNOSTIC USE

HALO® Patents: 7,354,729, 7,354,730 and other U.S. and international patents.

HemoGenix® and HALO® are registered trademarks of HemoGenix®, Inc

reagents produced by Promega Corporation and licensed under U.S. Patent HALO® was designed and developed by HemoGenix®, Inc and contains

Numbers 5,583,024, 5,674,713, 5,700,673 and others.

HemoGenix[®], Inc

No part of this instruction manual may be copied, duplicated or used without the express consent of HemoGenix®

For In Vitro Research Use Only. Not for clinical diagnostic use.

(Version 8-2013)

Manual

for Lympho-Hematopoietic Stem and Progenitor Cells

Suspension Expansion Culture[™] (SEC[™]) Master Mixes

CellExpand[™]

Web: www.hemogenix.com Fax: (719) 264-6253



Technical Support

Order online at www.hemogenix.com

Email: info@hemogenix.com

Fax: (719) 264-6253

Tel: (719) 264-6250

Ordering Information

Tel: (719) 264-6250

Email: info@hemogenix.com

1485 Garden of the Gods Road Colorado Springs, CO 80907 U.S.A. Tel: (719) 264-6250 HemoGenix[®], Inc Suite 152



TABLE OF CONTENTS



Step 2 - CellExpand [™] Cell Culture	 6. The CellExpandTM Protocol Step 1 - Cell Preparation A. Human, Non-Human Primate, Hose, Pig, Sheep or Canine Cells B. Rat or Murine Bone Marrow C. Isolation of Hematopoietic Subpopulations D. Cell Viability, Cell Counting and Cell Suspension Preparation E. Cell Thawing after Cryopreservation F. Mononuclear Cell Preparation 	5. Equipment, Supplies and Reagents Required, but not Provided	4. CellExpand [™] Reagents and Storage	3. Overview of the CellExpand TM Procedure	2. Intended Use: Applications and Advantages of Using CellExpand [™]	1. Introduction
11	ГГГГ 8 8 0 Q	6	С	4	2	

7. Recommendations and Tips Prior to Using CellExpandTM

12



.Ξ **Cell Suspension**

- any fractionation. a MNC preparation. However, murine bone marrow can usually be used without bone marrow and peripheral blood. Murine peripheral blood should also be human umbilical cord blood, human, primate, horse, pig, sheep, dog and rat The preferred cell suspension is a mononuclear cell suspension (MNC) of
- High concentrations of red blood cells can inhibit cell growth and expansion. It is therefore recommended to use an MNC preparation.
- are to be treated in cell culture, higher cell concentrations may be required. If cells have been treated (e.g. with cytotoxic drugs etc.) prior to cell culture, or

≣ Thawing and Storage of the CellExpand[™] Master Mix

- either in a 37°C water bath or at room temperature. Prior to using the CellExpand[™], remove the bottle from the freezer and thaw
- but can be refrozen. After use, CellExpand^m Master Mix is stable at 2-8°C for 1 month after thawing

1 **Humidity Chamber**

sufficient humidity to prevent evaporation. and drying out of the cultures. Even fully humidified incubators do not provide recommended to use an additional humidity chamber to prevent evaporation If microwell cultures or similar small volume cultures are performed, it is

1. INTRODUCTION

Changing the Paradigm

Assays You Can Trust Innovative Expertise You Can Count On

different stem cell populations and 7 progenitor cell populations from up to 8 species Suspension Expansion Culture[™] (SEC[™]) Master Mixes are available to culture up to 8 culture plates to large flasks, reactors and 3-dimensional culture systems. CellExpand™ progenitor cells in suspension culture using different culture vessels from microwell CellExpand[™] is used to culture and/or expand lympho-hematopoietic stem and

(human, non-human primate, horse, sheep, pig, dog, rat and mouse).

contains no growth factor cocktails, thereby allowing the user to add any growth at the required cell dose. cells to proliferate and expand. A CellExpand[™] Base Master Mix is also available which CellExpand[™] is a complete Master Mix. It contains everything required to allow the factors and/or cytokines at any concentration and for any species. Just add target cells



2. INTENDED USE: Applications and Advantages of Using CellExpand™

CellExpandTM is a complete reagent that can be used for stem cell research, basic and veterinary research and stem cell and regenerative cell therapy.

CellExpand[™] can be used for lympho-hematopoietic stem and progenitor cells from:

- Embryonic tissue
- Fetal tissue
- Spleen
- Bone marrow
- Peripheral blood
- Cord blood

derived from the following species:

- Human
- Non-human primate
- Horse
- Pig
- Sheep
- Sheep Dog
- Rat
- Mouse

CellExpand[™] is available for the following cell populations:

C = 11 T		Current Francisco Currentino a
Cell Type		טוטשנוו דמרנטו א/ כאנטאווופא
	Not specified	No growth factors or cytokines included
Stem Cells	L dS-ddH	IL-3, IL-6, SCF, TPO, Flt3-L
	2 dS-ddH	EPO, GM-CSF, IL2, IL-3, IL-6, IL-7, SCF, TPO, Flt3-L
	CFC-GEMM 1	EPO, GM-CSF, IL-3, IL-6, SCF, TPO, Flt3-L
	CFC-GEMM 2	EPO, GM-CSF, IL-3, IL-6, SCF, TPO
	CFC-GEMM 3	EPO, GM-CSF, G-CSF, IL-3, IL-6, SCF + TPO
	CFC-GEM 1	EPO, GM-CSF, IL-3, IL-6, SCF
	CFC-GEM 2	EPO, GM-CSF, IL-3, SCF
	CFC-GEM 3	EPO, GM-CSF, G-CSF, IL-3, SCF
Progenitor Cells	BFU-E 1	EPO, IL-3, SCF
	BFU-E 2	EPO
	GM-CFC 1	GM-CSF, IL-3, SCF
	GM-CFC 2	GM-CSF, G-CSF, IL-3, SCF
	GM-CFC 3	GM-CSF
	Mk-CFC 1	TPO, IL-3, SCF



treated with ACK lysis buffer.
 Removing excess density gradient reagent below the interface can result in

granulocyte and platelet contamination of the MNCs. In some cases, a trade-off of higher cell recovery with increased platelet contamination may be necessary. In these cases, cells can be centrifuged for 10 minutes at 400 x g instead of 300 x g, with no detrimental effect on the cells.

STEP 2. CellExpand[™] Cell Culture

Please refer to Section 7 for recommendations and tips prior to beginning this stage of the procedure.

Perform all cell culture procedures under sterile conditions in a biosafety cabinet.

With the exception of CellExpandTM that contains no growth factor cocktails, all CellExpandTM Master Mixes are complete and ready-to-use. No additional components should be added since these will dilute the components and result in non-optimized culture conditions.

The CellExpand[™] Method

- Transfer one or more bottles of frozen CellExpand[™] Master Mix to a 37°C water bath or thaw at room temperature.
- When thawed, mix the contents of bottle by inverting several times being careful not to cause bubbles.
- Prepare the cell suspension as required and adjust to the preferred working cell concentration.
- Prepare and label individual tubes or other culture vessels for each sample to be tested.
- 5. When the CellExpandTM Master Mix has thawed dispense the required amount of reagent into the culture vessel(s) being used <u>minus 10%</u>. **NOTE**: If using a CellExpandTM Assay Kit with no growth factors (Base Master Mix), first dispense the reagent followed by the specific growth factor(s) and or
- *cytokine(s).* Dispense the cell suspension at the correct working concentration so that 10% is added to the CellExpandTM Master Mix.

<u>о</u>

- Added to the CellExpand Master Mix. This is now the Culture Master Mix.
 Mix the Culture Master Mix thoroughly by vortexing and leave for a few minutes to settle.
- 8. Dispense the Culture Master Mix into the vessels
- Culture the cells at 37°C in a fully humidified incubator containing an atmosphere of 5% CO₂. If possible, use a 3-gas incubator to displace the atmospheric oxygen concentration (21%) to 5% O₂ with nitrogen. This increases the plating efficiency by reducing oxygen toxicity to the cells.
- 10. Lympho-hematopoietic stem cells initiate proliferation after about 24 hours and enter their exponential growth period after about 4-5 days. Progenitor cells usually demonstrate similar growth kinetics, but their total proliferation growth period will be shorter than for stem cells.

			 Specimens that have a high red cell concentrations can be separated again or
			additional cells.
			a sterile tube using a serological pipette and reserve until the cell count is
			If the cell separation is not clean and distinct, transfer the supernatant to
			spicules can be "picked out" using a manual 1ml pipette.
			in the cell interface. The fat will rise to the top during the second spin and the
			 Bone marrow often has fat globules and bone spicules; these can get caught
			with the supernatant after centrifugation.
	,		have some residual platelets; however, most of the platelets will be removed
nostic use or for cGMP conditions.	r clinical diagr	approved fo	sterile tube and centrifuge for 10 minutes @ 300 x g. The resulting cells will
arch use only and has not been	nd™ is for resea	CellExpar	 If fibring strands get caught up with the cells transfer entire supernatant to a
PORTANT:	MI		Do not discard the supernatant until a clean transfer has been made. Try to
			 Aspirate the supernatant using a serological pipette and save to a sterile tube.
c quantities.	10ml, 40ml and bulk	Available in	mononuclear cell laver. Care should be taken in this situation.
culture completion; (c) greater sensitivity.	tion; (b) more rapid	cell prolifera	 Eibrin strands can appear after separation. They are attached to the
actions providing (a) shorter lag times to initiate	allows for cell intera	 CellExpand[™] 	 Specimens that are not tresh or greater than 24 hours old, may not separate
essary. Sing using normal pipettes and liquid handlers	and accurate disnen	 Allows easy 	Density Gradient Centrifugation Troubleshooting
int species.	trom up to 8 differen	populations	
ent stem cell populations and 8 progenitor cell	detect up to 8 differ	Available to	improve viability.
no further additions. Just add cells.	aster Mixes require r	Pre-made M	Fractionating the MNCs by density gradient centrifugation will usually dramatically
		J	low proliferation ability.
M.	Jsing CellExpand ^T	Advantages of L	than 85%. Using cells with a viability lower than 85% will produce results with
	studies.	• Infindite cell	11. Perform a nucleated cell count and viability. The cell viability must be greater
dies.	of Action (IVIUA) stud	 Mechanism 	1ml manual pipette.
	lysis.	• Genetic anal	10. Add 1-2ml of IMDM and resuspend the cells, breaking up any clumps using a
	g studies.	Cell signaling	pellet.
	etry.	 Flow cytome 	9. Aspirate the supernatant after centrifugation taking care not to aspirate the cell
	on processes.	differentiatic	300 x g at room temperature.
ent screening for both proliferation and	and rapid high-cont	 Allows easy 	8. Add approx. 10-20 ml dPBS, mix gently and centrifuge the cells for 10 min at
ures.	3-dimensional cultu	reactors and	cells below the interface.
ulture vessel from microwell plates to cells	in virtually any cell c	 Cell growth i 	7. The vest of a particular the cells using a manual 1ml pinette. Do not remove
or cells in suspension culture.	f stem and progenit	 Expansion or 	7 Harvort the MNCs from the interface and transfer the cells to another starile
		h. h	layer above the MNC interface leaving approx. ¼" above the interface. Discard
	CellExpand™	Applications of (6. After centrifugation, remove the tube(s) gently and carefully aspirate the top
	tor. lL: Interleukin.	Colony Stimulating Fact	3. Cellfuluge for to min act, you x g of zo min action x g actionit temperature with NO brake.
proietin. SCF: Stem Cell Factor. TPO: Thrombopoietin. CSF:	rorming Cell. EPO: Erythro	B-lymphocyte Colony-Fi	allow the cells to mix with the reagent.
it - Erythroid. GM-CFC: Granulocyte-Macrophage Colony-	BFU-E: Burst-Forming Uni	Erythroid, Macrophage.	gradient reagent and cell suspension, the better the separation will be. Do not
stential - Stem and Progenitor. CFC-GEMM: Colony-Forming karyocyte. CFC-GEM: Colony-Forming Unit - Granulocyte,	P-SP: High Proliferative Po hroid, Macrophage, Megal	Abbrevitions used . HP Unit - Granulocyte, Erytl	a Pipette Aid on slow delivery. The cleaner the interface between the density
IL-7	B-CFC		 Using a sterile, serological pipette, dispense the diluted sample gently on top of the description reasons has believe the table of approx. All should be and the set of the
stimulator(s)			the density gradient reagent into a 50ml tube.
IL-2 (will usually require the addition of co-	T-CFC		aradient reagent into the tube. For samples greater than 3 ml. dispense 15ml of
Growth Factors/Cytokines	Cell Population	Cell Type	samples greater than 3ml, use a 50ml conical, plastic tube for separation.
Changing the Paradigm Innovative Expertise You Can Count O			Changing the Paradigm Innovative Expertise You Can Count On
MemoGenix® Assays You Can Trust			

10



3. OVERVIEW of the CellExpand[™] PROCEDURE

Using a CellExpand[™] is usually just a 2 step process

Step 1 – Cell Preparation

defined, pre-validated protocol to obtain a single cell suspension or procedures that are suggested in this manual. A dye exclusion viability and/or metabolic viability and Cells are not provided with CellExpandTM. Cells are prepared either by using a usernucleated cell count should be performed on all samples.

Step 2 – Cell Culture

cells should be transferred to a 37°C, humidified incubator gassed with 5% CO₂ and added at 10% of the final volume. Culture vessels containing lympho-hematopoietic specific lympho-hematopoietic cell populations. The reagent is usually provided in preferably 5% O_2 . 10ml or 40ml bottles. Depending on the volume of CellExpand[™] to be used, cells are CellExpand[™] is a complete, ready-to-use Master Mix for growing and/or expanding

E. Cell Thawing after Cryopreservation

Thawing of Cells

or alleviate the possibility of clumping during cell thawing, it is recommended that was a MNC or similar fraction, the chances of clumping will be low. However to reduce and other cell components will rupture and release DNA. Large amounts of released cells, granulocytes and other cell populations. When the cells are thawed, granulocytes aliquots of cells only (1-1.5ml). DNA will clump together encasing cells. If the cell preparation originally cryopreservec DNase be added to the cell suspension. The following procedure is used for small If cells are cryopreserved as a total nucleated cell product, they will contain red blood

- Thawing of the vial contents is initiated in a 37°C water bath, by swirling the vial for approx. 1 min.
- Ņ and carefully unscrew the vial lid. the water bath, sterilized the outside of the vial by spraying with 70% ethanol When a small ball of ice still remains in the vial (1-2 min), remove the vial from
- ω to the total volume in the vial to achieve a concentration of 6µg/ml before It is possible that clumping can occur at this stage, in which case, add DNase proceeding to the next step.
- 4 stage. In this case, DNase at a final concentration of 6µg/ml should be added added to this 20ml of thaw medium. However, clumping can also occur at this Using a 1ml pipette, gently mix the contents of the vial and transfer to a 50ml tube containing 20ml of thaw medium. Up to 3 vials of the same cells can be before proceeding to the next step.
- Ś pipetting to mix the cells. This could cause further rupture of cells and the Gently mix the cells by swirling the contents of the tube. Do not use repeat release of DNA resulting in increased clumping.
- Centrifuge the cells at 300 x g for 10 min at room temperature and discard the supernatant after centrifugation.

<u>о</u>

2 Resuspend the cells in 1ml of culture medium. If necessary, add 6µg/ml DNase.

F. Mononuclear Cell Preparation

Density Gradient Centrifugation

centrifugation. It should be emphasized that when the red blood cell content is high, cases, however, density gradient centrifugation will produce a high quality MNC even density gradient centrifugation may not provide optimal depletion. In most Red blood cells and other contaminating cells can be removed by density gradient fraction and will also increase the viability of the cell sample.

density gradient centrifugation using NycoPrep is as follows: centrifugation reagent is brought to room temperature prior to use. The procedure for to cultures without any deleterious effects. It is essential that the density gradient procedure. This is because all Ficoll-derived reagents are toxic to cells. It is HemoGenix[®] does not recommend using a Ficoll density gradient centrifugation recommended to use NycoPrep 1.077, which is non-toxic and can also be added

- 2 Dilute the sample to be fractionated with an equal volume of Dulbecco's phosphate buffered saline (dPBS) and mix gently by inversion
- For samples 3ml or less use a 15ml conical plastic tube for separation. For

 Assays You Can Tust Changing the Paradigm Innovative Expertise You Can Count On tube. Flush out the marrow through the bone cavity and withdraw part of the cell 	4. CellExpand TM Reagents and Storage
suspension through the bone and into the syringe. 7. Flush the cell suspension through the bone and repeat steps 6 and 7 two to three times. When finished, the bone should appear translucent, indicating that most of the cells have been flushed out of the cavity.	CellExpand [™] is usually supplied as 10ml or 40ml volumes of complete, use Master Mix containing all the components necessary to grow and/v the required cell population(s). The reagent is shipped complete and n
 Remove the empty bone and replace it with the next bone until the marrow from all bones has been flushed out of the cavities. Let the cell suspension settle for 1-2 minutes to allow large debris to fall to the 	CellExpand [™] reagents are shipped frozen. Upon arrival, transfer the rea
bottom of the tube. 10. Using a small gauge (22-25) needle and syringe, slowly withdraw the cell suspension leaving the large debris in the tube and transfer it to a new tube,	וופבצבו טוונו טאבט. Once thawed, CellExpand™ can be stored at 4°-8°C for 1 month, but cal rafrozen and thawed at a later time
11. If necessary, add medium to achieve the required volume.	Growth factors and/or cytokines should only be added if the No Grow
To further purify a cell suspension, it is recommended to underlay 1-2ml of NycoPrep 1.077 (Axis-Shield) under the cell suspension in a sterile, 15ml conical tube with screw can and centrifuce the cells according to the manufacturer's protocol	CellExpand [™] Base Master Mix was ordered.
C. Isolation of Hematopoietic Subpopulations Providing sufficient cells are available, subpopulations of stem and progenitor cells can be isolated and purified for use with CellExpand [™] . HemoGenix [®] recommends using magnetic cell isolation procedures (e.g. Miltenyi Biotech), since these allow rapid isolation of different cell populations with substantial purity, viability and yield. Please see table of recommended cell concentrations to determine the optimal, final cell concentration to use with CellExpand [™] reagents.	Good laboratory practices and universal protective precautions s undertaken at all times when handling the kit components as we tissues. Material safety data sheets (MSDS) are included in each li
 D. Cell Viability, Cell Counting and Cell Culture Suspension Preparation 1. For dye exclusion viability methods, use trypan blue and a hemocytometer or automated method, flow cytometry using 7-AAD or another vital stain. Note that dye exclusion viability methods detect membrane integrity. They do not detect cellular and mitochondrial integrity and therefore metabolic viability. A viability of 85% or greater should be obtained when using dye exclusion viability methods only. It is recommended not to use cell suspensions with a viability of less than 85% since these cells will not be able to sustain proliferation 	
 Determine the cell concentration using either a hemocytometer or electronic cell/particle counter. Adjust the cell suspension concentration to that recommended in Table 1. Note the working cell concentration per ml is 100 x the final cell concentration per 	
 may be required. Prepare the total volume of cell suspension required using IMDM or PBS. The volume of the adjusted cell suspension required will be 10% of the total volume of the adjusted cell suspension required will be 10% of the total volume of the CellExpand[™] Culture Master Mix. 	



5. EQUIPMENT, SUPPLIES AND REAGENTS REQUIRED, BUT NOT PROVIDED

The following is a selection of equipment, supplies and reagents that might be required for use with CellExpandTM.

Equipment and Supplies

- 1. 35mm Petri dishes that are not tissue culture treated.
- 2. Inverted microscope fitted with 10x oculars and a 4x and 10x objectives
- 3. Laminar Flow Biohood.
- 4. Sterile plastic tubes (5ml, 10ml).
- 5. Single channel pipettes, preferably electronic (e.g. ViaFlow or Rainin EDP pipettes for variable volumes between 1µl and 1000µl).
- Sterile pipette tips.
 Repeater pipette w
- Repeater pipette with positive displacement (e.g. Eppendorf Repeater Strean, Oxford, Gilson Distriman or Rainin AutoRep E) for dispensing CellExpand^m reagent.
- 8. Sterile syringes for repeater pipette
- 9. Vortex mixer.
- Tissue culture incubator, humidified at 37°C with 5% CO₂ (minimum requirement) and 5% O₂ (preferable).
- Hemocytometer or electronic cell counter to determine cell concentration
 Flow cytometer or hemocytometer for determining viability.

Reagents

- 1. Sterile Phosphate Buffered Saline (PBS)
- Iscove's Modified Dulbecco's Medium (IMDM)
- 3. Density-gradient medium (e.g. NycoPrep 1.077, Axis-Shield).
- 7-AAD, propidium iodide or trypan blue for viability assay.

4

6. The CellExpand[™] PROTOCOL

Changing the Paradigm

Assays You Can Trust Innovative Expertise You Can Count On

PLEASE READ THE FOLLOWING PROTOCOL CAREFULLY. SEE SECTION 7 BEFORE PERFORMING THE ASSAY

Good laboratory practices and universal protective precautions should be undertaken at all times

Using CellExpand[™] is usually a 2 step process.

- Step 1 Cell preparation.
- Step 2 − CellExpand[™] culture setup

All steps must be performed in a laminar flow biohazard hood

STEP 1 – Cell Preparation

A. Human, Non-Human Primate, Horse, Pig, Sheep or Canine Cells

- 1. For best results, CellExpand[™] requires that target cells be separated from red blood cells since these can interfere with cell growth and expansion. Neutrophils and platelets should also be removed. Therefore, the starting cell suspension should be a mononuclear cell (MNC) suspension. HemoGenix[®] recommends separating the MNCs using NycoPrep 1.077 (Axis-Shield) density gradient centrifugation media. This separation procedure should be used for human, non-human primate, horse, pig, sheep and canine cells. Ficoll-Paque can also be used, although this is toxic to cells. Follow the manufacturer's protocol to prepare the MNCs.
- 2. Resuspend the cells in IMDM or PBS

Human umbilical cord blood.

It is recommended to deplete human umbilical cord blood of erythrocytes using a current Hetastarch[®] protocol or a density gradient separation per the manufacture's protocol so that the majority of the erythrocytes are removed from the cell suspension. It is not recommended to perform the CFC assay with a high concentration of erythrocytes present for the same reasons as above.

B. Rat or Murine Bone Marrow

- 1. Remove organs (femora and tibia (optional)) under aseptic conditions.
- Remove as much muscle from the bones as possible.
- Using a sterile blade, first cut off the proximal (hip joint) end below the ball joint at right angles to the longitudinal length of the bone. Then cut off the distal end (above the patella or knee).
- Transfer sufficient sterile medium to a tube so that it will cover the whole bone, approximately 1-2ml. (Some of the medium provided with the kit can be used for this purpose).
- 5. Half fill a syringe (1-3ml) with sterile medium and, using a needle gauge that will enter the bone cavity without cracking the bone, insert the needle into the proximal end and immerse the whole bone in the medium contained in the